



Natural Tooth Color Ceramic Feldspathic Veneer For Anterior And Posterior Teeth

Our Product Introduction

for more products please visit us on adldentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Esthetics: Highly Aesthetic And Natural Appearance
- Durability: Long-lasting
- Customization: Customizable Shape And Size
- Feature: Non-invasive And Minimal Prep Is Required
- Surface Texture: Smooth
- Color: Natural Tooth Color
- Application: Suitable For Anterior And Posterior Teeth
- Material: Feldspathic Ceramic
- Highlight: **Ceramic Feldspathic Veneer, Natural Tooth Color Feldspathic Veneer, Posterior Teeth Feldspathic Veneer**



More Images



Product Description

Smooth Surface Textured Natural Tooth Color Feldspathic Ceramic Veneer

Product Description:

One of the standout features of Feldspathic Veneers is their long-lasting durability. They are designed to withstand the daily wear and tear of normal use, making them a great investment for anyone looking to improve their smile.

In addition to their durability, Feldspathic Veneers are known for their highly aesthetic and natural appearance. They are carefully crafted to blend in seamlessly with the surrounding teeth, giving patients a beautiful and natural-looking smile.

Another benefit of Feldspathic Veneers is that they are non-invasive, with minimal prep required. This means that patients can enjoy the benefits of a beautiful smile without having to undergo extensive dental work.

Overall, Feldspathic Veneers are an excellent choice for anyone looking to improve their smile. With their natural appearance, long-lasting durability, and minimal prep requirements, they are a great investment for anyone looking to achieve a beautiful, healthy-looking smile.

Features:

Product Name: Feldspathic Veneer

Color: Natural Tooth Color

Durability: Long-lasting

Translucency: High Translucency

Thickness: 0.3-0.7mm

Material: Feldspathic Ceramic

Our Feldspathic Veneer is a type of Porcelain Veneer made from Feldspathic Ceramic. It provides a natural tooth color and high translucency, while also being long-lasting. With a thickness of 0.3-0.7mm, it is a great option for anyone looking for Feldspathic Porcelain Veneers or Porcelain Veneers in general.

Technical Parameters:

Product Attribute	Description
Product Name	Feldspathic Porcelain Veneers
Thickness	0.3-0.7mm
Customization	Customizable Shape And Size
Durability	Long-lasting
Surface Texture	Smooth
Material	Feldspathic Ceramic
Application	Suitable For Anterior And Posterior Teeth
Esthetics	Highly Aesthetic And Natural Appearance
Color	Natural Tooth Color
Feature	Non-invasive And Minimal Prep Is Required
Translucency	High Translucency

Applications:

These veneers are suitable for both anterior and posterior teeth, making them ideal for a wide range of dental applications. The smooth surface texture of the veneers ensures a natural and seamless look, and their thickness can range from 0.3mm to 0.7mm depending on the patient's needs.

ADL Feldspathic Veneers are designed to be long-lasting, meaning that they will provide patients with a beautiful and functional smile for years to come. With a supply of 630 technicians, these veneers can be customized to fit any patient's unique needs and preferences. Overall, ADL Feldspathic Veneers are a high-quality and reliable choice for cosmetic dentistry. With their durability, customization options, and smooth surface texture, they are sure to provide patients with a beautiful and natural-looking smile. Payment terms are flexible, with both Paypal and T/T options available.

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong
Province, P.R.China