




Titanium Sandblasted Implant Locator Bar Silver Gold Dental Implant Locator Attachment

Our Product Introduction

for more products please visit us on adidentalab.com

Basic Information

- Place of Origin: Shenzhen, China
- Brand Name: ADL
- Certification: ISO/CE/FDA
- Minimum Order Quantity: Negotiable
- Price: Negotiable
- Packaging Details: Plastic bags and box
- Delivery Time: 72 Hours
- Payment Terms: Paypal or T/T
- Supply Ability: 630 Technicians



Product Specification

- Connection Type: Internal Hex
- Color: Silver And Gold
- Attachment Type: Locator
- Number Of Attachments: 2 Or More
- Sterilization Method: Autoclave
- Surface Treatment: Sandblasting
- Material: Titanium
- Compatibility: Implant Supported Overdentures
- Highlight: **adl Titanium Sandblasted Implant Locator Bar, adl Sandblasted Implant Locator Bar, Titanium Sandblasted Implant Locator Bar Silver**



More Images



Product Description

Titanium Sandblasted Locator Bar Attachments in Silver and Gold Implant bar and Locator Attachment

Product Description:

The LOCATOR with milled implant bar is an excellent choice for implant-supported overdentures. The attachment type is a Locator, which allows for easy insertion and removal of the denture for cleaning and maintenance. The product is made of high-quality titanium, which is known for its strength, durability, and biocompatibility, making it an ideal material for dental implants.

The LOCATOR with milled implant bar provides exceptional stability and support for the overdenture. This combination of the LOCATOR Attachment and milled implant bar creates a secure and stable platform for the denture, which can help to improve patient comfort and confidence while also reducing the risk of implant failure.

Features:

Product Name: Implant bar and Locator Attachment

Attachment Type: Locator

Connection Type: Internal Hex

Color: Silver And Gold

Compatibility: Implant Supported Overdentures

Surface Treatment: Sandblasting

Some possible variations of the product name include:

LOCATOR Implant Attachment

LOCATOR with milled implant bar

Technical Parameters:

| Product Attribute | Value |
|-----------------------|------------------------------------|
| Product Name | Implant bar and Locator Attachment |
| Number Of Attachments | 2 Or More |
| Compatibility | Implant Supported Overdentures |
| Sterilization Method | Autoclave |
| Color | Silver And Gold |
| Connection Type | Internal Hex |
| Attachment Type | Locator |
| Surface Treatment | Sandblasting |
| Material | Titanium |

Applications:

The Implant bar with Locator Attachment is perfect for use in different dental applications. It is ideal for patients who require a reliable and stable fitting for their dentures. This product is suitable for patients with partial or full dentures, ensuring that they are comfortable and can chew and speak with ease.

The LOCATOR Implant Attachment is also perfect for patients who have undergone a bone graft or have limited bone density. It provides a secure fit for the dentures, ensuring that the patient experiences minimal movement of the dentures. This reduces the likelihood of discomfort and gum irritation, which is common with traditional dentures.

The Implant bar with Locator Attachment is delivered in plastic bags and boxes, ensuring that it is safe during transportation. The product has a fast delivery time of 72 hours, making it perfect for dental professionals who want to provide their patients with a quick solution for their dental needs.

The Implant bar and Locator Attachment is available at a negotiable price, and the minimum order quantity is also negotiable. The product can be purchased using PayPal or T/T, and the supply ability is 630 technicians, ensuring that dental professionals can order the product in bulk.

In conclusion, the Implant bar with Locator Attachment is a reliable and long-lasting solution for dental professionals and their patients.

The product is made from high-quality titanium material, making it durable and long-lasting. The LOCATOR Implant Attachment is perfect for patients who require a secure and comfortable fit for their dentures, ensuring that they can chew and speak with ease. The product is available at a negotiable price, making it affordable for dental professionals, and has a fast delivery time of 72 hours.



Advanced Dental Laboratory

+86 18665847934

info@adidentalab.com

adidentalab.com

Plant 1, Fuqiao District 2, Qiaotou Community, Fuhai Street, Baoan District, Shenzhen City, Guangdong Province, P.R.China